

panmed | **u.s.**
PAN MEDICAL | US CORPORATION

RADIOPAQUE BONE CEMENT FOR
VERTEBROPLASTY AND KYPHOPLASTY

TECHNICAL REPORT

DESCRIPTION

OSTEOFLEX® is a PMMA bone cement, specially designed for percutaneous bone augmentation procedures.

Its composition has been optimized to give the best properties to the surgeons to treat patients in a safety way.

After a fast mixing time, **OSTEOFLEX®** is ready to use in a paste form.

It is a medium viscosity bone cement with a very low exothermicity.

The graphic features a light blue background with a network of darker blue, branching lines resembling a biological or cellular structure. The text 'OsteoFlex' is prominently displayed in a large, bold, dark blue font, with 'Radiopaque bone cement' written below it in a smaller, dark blue font.

OsteoFlex
Radiopaque bone cement



INDICATIONS

OSTEOFLEX® is indicated for the fixation of pathological fractures of the vertebral body using vertebroplasty or kyphoplasty procedures.



- Painful osteoporotic vertebral compression fractures
- Benign symptomatic lesions (haemangioma)
- Malignant lesions (tumors, metastasis lesions, myeloma)



CLINICAL
DATA



Treatment of vertebral myeloma



Treatment of L2 vertebral myeloma




Viscosity

The viscosity curve of **OSTEOFLEX®** follows a particular profile with a long working time up to 20 minutes.

The viscosity kinetics and setting time of **OSTEOFLEX®** have been specifically designed for a safe application.

It is a medium viscosity bone cement able to suitably flow from the injection system through a trocar into the bone.

This special and stable viscosity over a period of 20 minutes is largely enough to perform a continuous radiological injection checking.




*Viscosity profile of **OSTEOFLEX®** at 20°C*



APPLICATION TIME

The long working time of **OSTEOFLEX®** (20 minutes) is sufficient to:


- Inject the cement into the vertebral body
- Take time to check via different views
- Work comfortably under CT scan
- Perform a multi level injection procedure.





OSTEOFLEX® has a high rate of zirconium dioxide opacifier for an ideal and maximum visibility during injection.

It can be used either under fluoroscopy or CT scan.




OSTEOFLEX® radiopacity

LOW EXOTHERMIC REACTION

The maximum temperature during the setting is about 45°C (according to ISO5833 standard).

This extraordinary low exothermic reaction of **OSTEOFLEX®** prevents the risk of thermic necrosis of bone tissue, in comparison with standard bone cements.



PREPARATION AND INJECTION

Mixing phase


- The two components (powder and liquid) have to be mixed correctly for 30 seconds in the shaker.

Filling the syringe

- Transfer the cement into syringes.

Injecting phase (up to 20 minutes)

- No waiting time is required. Immediately after filling the injector (approx. 5 minutes after beginning the mixing process), the cement becomes a paste, ready to be injected safely.
- This paste constitution corresponds to the best viscosity preventing leakage.
- Still injectable up to 20 minutes at 20°C.
- The injection of cement must take place under continuous radiological guidance. When the operator estimates that the vertebral filling is satisfactory or when a risk of leak of cement becomes apparent, the injection can be stopped.



- Mixing
- Injection time/working time *
- Injector filling and waiting time
- In situ hardening time **

As usually recommended, when the cement injection is finished, keep the patient in place for 15 minutes to get a complete in situ polymerization of the cement.




**PRODUCT
CLASSIFICATION**

Class: IIb medical device

Applicable directive: 93/42/EEC

OSTEOFLEX® comes in the form of a sterile liquid ampoule of monomer and a sterile powder pouch of powder polymer.



MEDICAL DEVICES

PACKAGING AND STERILISATION

The liquid in the ampoule is sterilized by ultra-filtration, and the ampoule blister is sterilized using ethylene oxide. The powder is inside a double sterile pouch. This double pouch is sterilized by 25 kGy gamma rays.

Disposable, single use only.






Powder

Poly methyl methacrylate and copolymer
Benzoylperoxide
Zirconium Dioxide

Liquid

Methyl methacrylate
N-N dimethyl-p-toluidine
Hydroquinone



Pan Medical Ltd
Barnett Way
Barnwood
Gloucester
GL4 3RT

www.pan-medical.co.uk

panmed | **u.s.**
PAN MEDICAL | US CORPORATION